




*medical*

---

## **Blood and Fluid Warming Inserts**

### **Compatible for use with most Dry Heat Fluid Warmers**

- Optional 3 way tap and extension available
- Guaranteed free of Latex
- Extended patient delivery line for ease of access during procedures
- Additional clamp on inlet line to seal the line before disposal
- Locking Luer fittings on both inlet and outlet lines for enhanced security
- Location holes punched through to allow easy fitment to the warmer
- Clear and rigid bubble trap to allow easy visual check of fluid flow with minimal risk of chamber collapse during extended procedures
- CE marked under G.M.P. Certification No. CE 0482 to comply with the latest regulations for medical devices. Council Directive No. 93/42/EEC
- Prices dependant on quantity and/or annual useage
- Standard insert is code No. Tinsert/Std (T105100) and the 3 way tap insert is code No. Tinsert/3WT (T105100SC) both come in cartons of 4 x 10 pieces



---

## **JMW Medical Ltd**

Pentland Industrial Estate, Loanhead, Midlothian, EH20 9QH

Tel: 0131 448 0205 Fax: 0131 448 0208

Web: <http://medical.jmw-group.co.uk/> email: [medical@jmw-group.co.uk](mailto:medical@jmw-group.co.uk)